

## PRODUCT BROCHURE

---

# KERN KUP-06

Date: 27.08.2025

ENGLISH



# Image not available

Interface adapter with cable. Interface with Bluetooth BLE (v4.0)

KUP, cable length 0,15 m

---

### CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- The KERN Universal Port (KUP) enables the exchange of data and control commands between balances and the network or PC, mobile end unit, printer, etc. using several interfaces: RS-232, USB device, Ethernet, Bluetooth and WiFi. In this way we ensure the perfect integration of your KERN balance into production or processing chains, which gives the benefit of a simple, consistent workflow
- In accordance with the „Plug and Play“ principle, the KUP for RS-232 or USB, Ethernet, Bluetooth or WiFi will only be inserted into the DB 15 socket which is integrated into the device. It is therefore possible to retrofit or replace interfaces conveniently without having to open the balance casing and no complicated installation. This option therefore offers a high level of flexibility at low additional cost

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	KUP-06
Model-Serie	KUP
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	4045761253957
Customs tariff number	90160090

Packaging component - 27 g  
by weight - styrofoam

Delivery time 1 d

### FUNCTIONS

Standard

### Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	225×145×40 mm
Gross weight	0.15 kg
Net weight	0.10 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	27 g