

PRODUCT BROCHURE

**KERN
3 17-07-2022i**

Date: 21.07.2025

ENGLISH**Image not
available**

Individual weight, OIML Class E2, 100 g, Knob, stainless steel polished

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com

Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	317-07-2022i
Model-Serie	317
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Trade parts list element
EAN code (Model)	4045761245334
Predecessor Model	317-07; 317-07-2021a
Successor Model	317-07
Customs tariff number	84239090
IsKERNReferenceArticle NonMandatory	False

Technical data - Primary

Nominal value	100 g
Tolerance (OIML)	0,16 mg
Case aluminium	317-070-600
Case plastic	317-070-400
Case wood	317-070-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

OIML class	E2
OIML conform	yes
Density (kg/m ³)	7950

Volume (cm³) 12,5

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	22×22×38 mm
Gross weight	0.10 kg
Net weight	0.10 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-337
-------------------	---------

FUNCTIONS

Standard

ENGLISH

ACCESSORIES

Model	Description
317-070-600	Aluminium Protective Box KERN 317-070-600
317-070-100	Wooden Box for Set of Weights KERN 317-070-100
317-070-400	Plastic Weight Box KERN 317-070-400

SERVICES

Model	Description
921-337	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-337
962-337	DAkkS Calibration KERN 962-337
952-337	Verification KERN 952-337
962-337R	DAkkS Recalibration KERN 962-337R