

**PRODUCT BROCHURE**

---

**KERN  
KUM-01**

Date: 26.08.2025

ENGLISH



**Image not  
available**

Internal data interface RS-232

---

**CONTACT**

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

Phone : +49 7433 9933-0  
Fax : +49 7433 9933-149  
Email : [info@kern-sohn.com](mailto:info@kern-sohn.com)  
Web : [www.kern-sohn.com](http://www.kern-sohn.com)



Scan here for more  
information

## DESCRIPTION

- The KERN Universal Modules (KUM) enable the exchange of data and control commands between balance and network or PC, mobile end unit, printer etc. using several interfaces: RS-232, USB device, Ethernet, Bluetooth and WiFi. In this way we ensure the perfect integration of your KERN balance into production or processing chains, which gives the benefit of a simple, consistent workflow
- The KUM interface modules for RS-232, USB, Ethernet, Bluetooth and WiFi are plugged straight onto the board to suit individual choice and this means that the IP protection of the housing remains guaranteed

## SPECIFICATIONS

### Article codes, product groups & Status information

Model article-code	KUM-01
Model-Serie	KUM
Assortment	KERN
Product Type	General weighing accessories
Product group	Accessories (printer, etc.)
Product usage type	Accessories/DL
EAN code (Model)	404576 1476738
Customs tariff number	90160090

Delivery time 1 d

## FUNCTIONS

Standard

## Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	225×145×40 mm
Gross weight	0.50 kg
Net weight	0.30 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	100 g
Packaging component - by weight - styrofoam	100 g